

Houndstooth Bolster Furniture Protector



Using Your Bolster Furniture Protector & Washing Instructions

Thank you for purchasing the Bolster Furniture Protector. Your four legged best friends will enjoy sleeping on the couch, and your couch will remain clean and protected for years to come. The protector is fast and easy to put on, take off, and machine wash and dry.

Step 1:

Remove the protector cover and filled pillow from the plastic bag. Confirm that you have received 1) a protector cover (the bolster is empty), and 2) a filled white bolster pillow.

Step 2:

Open the bolster zipper on the protector cover. Put the filled white bolster pillow into the bolster. Make sure you push the filled pillow right to the end of the bolster. Close the zipper.

Step 3:


Locate the rubber "handles" on the underside of the protector. These handles help keep the protector from moving on the couch. Push the handles down between the seat cushions and back rest.

Step 4:

This is the most important step - enjoy your four legged best friend sleeping comfortably on your couch (or is it 'their' couch?) knowing your furniture will be protected.



Rubber Handle
This is what you push down between the cushions to make sure the protector stays in place.



Washing Instructions:

Your Furniture Protector is fast and easy to clean.

Wash and dry separately from other items. The filled pillow within the bolster must be removed prior to washing. Do not wash the pillow.

Wash the protector with cold water on a gentle setting using a mild laundry detergent, such as Woolite.

Tumble machine dry on a cool to warm setting.

Do not use bleach, and do not dry clean as these will damage the upholstery grade fabric. In the event of a pet 'accident' remove the protector immediately from the couch. We recommend using an enzyme based urine treatment directly on the accident. Follow the directions on the urine treatment product.

